

Gallium-67 Radiochemical Gallium Chloride Solution

Gallium (Ga-67) Chloride Half-Life: 3.26 days

GALLIUM-67 RADIOCHEMICAL GALLIUM CHLORIDE SOLUTION

Gallium chloride is usually complexed with sodium citrate to produce gallium citrate. It is useful for imaging various inflammatory conditions and tumours such as Hodgkin's disease, lymphoma (with the exception of lymphocytic), hepatoma and bronchogenic carcinoma.

Nordion has an unparalleled base of technical expertise and production capabilities. State-of-the-art facilities provide a reliable supply of isotopes to customers around the world. Our strong focus on customer service and quality has made us a world leader.

PROCESS DATA

Nuclear Reaction

$^{68}\text{Zn} (p,2n) ^{67}\text{Ga}$

Assay

High resolution gamma spectrometry at 184.6 keV

PRODUCT SPECIFICATION

Chemical Form

Gallium chloride in dilute HCl

Normality

0.1 N \pm 0.05 N HCl

Specific Activity*

≥ 1000 Ci/g

Activity Concentration*

≥ 40 mCi/mL

Radiopurity*

$^{67}\text{Ga} \geq 99.4\%$

$^{66}\text{Ga} \leq 0.5\%$

Others $\leq 0.1\%$

Chemical Purity*

Fe, Zn, Cu each ≤ 0.2 $\mu\text{g}/\text{mCi}$ ^{67}Ga

*Activity and radiopurity specifications are at noon, Pacific Time, one day after shipment.

Please note that the radiochemicals listed are not tested for sterility or pyrogenicity. Verification of their suitability for use in humans is the sole responsibility of the purchaser.

PACKAGING

- Hypovial – crimp cap glass vial with septum (6 mL, 10 mL)
- Microvial (1 mL, 3 mL)
(Other vial options are available).

AVAILABILITY

Dispensed daily

PRODUCT ORIGIN

Vancouver, Canada

WORLD CLASS DELIVERY SERVICE

Nordion has a world class distribution system that includes a unique, software-based order management system. Carefully managed transportation logistics ensure products reach their destination on time.



Gallium-67 Decay Table Half-Life: 3.26 days

HOURS	0	1	2	3	4	5	6	7	8	9	10	11
0	1.0000	0.9912	0.9824	0.9738	0.9652	0.9567	0.9482	0.9399	0.9316	0.9234	0.9152	0.9071
0.5	0.8991	0.8912	0.8834	0.8756	0.8678	0.8602	0.8526	0.8451	0.8376	0.8302	0.8229	0.8157
1	0.8085	0.8013	0.7943	0.7873	0.7803	0.7734	0.7666	0.7598	0.7531	0.7465	0.7399	0.7334
1.5	0.7269	0.7205	0.7142	0.7079	0.7016	0.6954	0.6893	0.6832	0.6772	0.6712	0.6653	0.6594
2	0.6536	0.6478	0.6421	0.6365	0.6309	0.6253	0.6198	0.6143	0.6089	0.6035	0.5982	0.5929
2.5	0.5877	0.5825	0.5774	0.5723	0.5672	0.5622	0.5573	0.5524	0.5475	0.5427	0.5379	0.5331
3	0.5284	0.5238	0.5191	0.5146	0.5100	0.5055	0.5011	0.4966	0.4923	0.4879	0.4836	0.4794
3.5	0.4751	0.4709	0.4668	0.4627	0.4586	0.4545	0.4505	0.4466	0.4426	0.4387	0.4348	0.4310
4	0.4272	0.4234	0.4197	0.4160	0.4123	0.4087	0.4051	0.4015	0.3980	0.3945	0.3910	0.3875
4.5	0.3841	0.3807	0.3774	0.3740	0.3707	0.3675	0.3642	0.3610	0.3578	0.3547	0.3516	0.3485
5	0.3454	0.3423	0.3393	0.3363	0.3334	0.3304	0.3275	0.3246	0.3217	0.3189	0.3161	0.3133
5.5	0.3105	0.3078	0.3051	0.3024	0.2997	0.2971	0.2945	0.2919	0.2893	0.2867	0.2842	0.2817
6	0.2792	0.2768	0.2743	0.2719	0.2695	0.2671	0.2648	0.2624	0.2601	0.2578	0.2556	0.2533
6.5	0.2511	0.2488	0.2467	0.2445	0.2423	0.2402	0.2381	0.2360	0.2339	0.2318	0.2298	0.2278
7	0.2257	0.2238	0.2218	0.2198	0.2179	0.2160	0.2141	0.2122	0.2103	0.2084	0.2066	0.2048
7.5	0.2030	0.2012	0.1994	0.1977	0.1959	0.1942	0.1925	0.1908	0.1891	0.1874	0.1858	0.1841
8	0.1825	0.1809	0.1793	0.1777	0.1762	0.1746	0.1731	0.1715	0.1700	0.1685	0.1670	0.1656
8.5	0.1641	0.1627	0.1612	0.1598	0.1584	0.1570	0.1556	0.1542	0.1529	0.1515	0.1502	0.1489
9	0.1475	0.1462	0.1450	0.1437	0.1424	0.1412	0.1399	0.1387	0.1375	0.1362	0.1350	0.1338

The half-life values stated are as of time of publication. Nordion values may vary.

CORPORATE HEADQUARTERS

447 March Road
 Ottawa, ON, Canada K2K 1X8
 Tel: +1 613 592 2790
 Fax: +1 613 592 6937

REGIONAL OFFICE

4004 Wesbrook Mall
 Vancouver, BC, Canada V6T 2A3
 Tel: +1 604 228 8952
 Fax: +1 604 228 5990

ASIA PACIFIC SALES OFFICES

Hong Kong

9/F, 22 Yee Wo Street
 Causeway Bay, Hong Kong
 Tel: +852 2827 8666
 Fax: +852 2827 8302

Japan

Room 905, Tokyo Royal Plaza
 1-18-11, Uchi-kanda, Chiyoda-ku
 Tokyo 101-0047, Japan
 Tel: +81 3 5283 6872
 Fax: +81 3 5283 6873

www.nordion.com